

Volcano R-100™

Reusable Pullback Device



Precise and Reusable

- ▶ Calibrated to ensure a smooth, accurate pullback rate: 0.5 mm/sec or 1.0 mm/sec
- ▶ Fully reusable device
- ▶ Adjustable, secure catheter grip
- ▶ Custom sterile drape available for quick and easy setup

Intended to be used with the Volcano family of digital intravascular imaging catheters

R-100 Technical Specifications

Height	Width	Depth	Weight	Power
9.9 cm (3.9")	20.3 cm (8")	2.5 cm (1")	0.95 kg (2.1 lb)	9V Alkaline, replaceable

Ordering Information

- ▶ Catalog number: R100
- ▶ Custom Drapes: PBKBG-10



PRECISION GUIDED THERAPY

Volcano Corporation

Corporate Headquarters
3721 Valley Centre Drive, Suite 500
San Diego, CA 92130 USA
Phone: 800-228-4728
Fax: 858-720-0325

Volcano Europe BVBA/SPRL

EU Authorized Representative
Excelsiorlaan 41
B-1930 Zaventem Belgium
Phone: +32 2 679 1076
Fax: +32 2 679 1079

Volcano Japan

Japan Headquarters
Hamamatsucho Square 6F
1-30-5, Hamamatsucho,
Minato-ku, Tokyo 105-0013, Japan
Phone: +81-3-6430-9400
Fax: +81-3-6430-9401

For further information about Volcano and its products, please visit www.volcanocorp.com.
Volcano and the Volcano logo are registered trademarks of Volcano Corporation.
R-100 is a trademark of Volcano Corporation.

601-0500.07/003